



## Your first step to become a nymph

Your path to become beautiful like a nymph starts here, with **CLIZIA RF**, the cutting-edge system for multipolar radiofrequency treatments for face and body. Once you have prepared the skin and chosen the cosmetic, you are ready to start the sweetest beauty experience. You can customize the session according to the clients needs. Everything can be adjusted: program, handle, timing, intensity, always according to the skin area and the skin sensitivity. For custom-cut treatments. Advanced technology (Multipolar RF) and power control make CLIZIA RF safe and effective from the first sessions: face and body appear immediately more tonic, firm and compact and the skin complexion brighter. All ensuring the best comfort. CLIZIA is not only skilful but it is also beautiful: stylish design and soft lines make the device sensual like Clizia, the charming nymph turned into sunflower in love with the god of sun Apollo. And thanks to the high performance of CLIZIA your client will change shape as well, blossoming again in all her beauty.

### How it works

CLIZIA offers firming treatments for face and body by means of **radiofrequency at 0,5 MHz e 1 MHz**, a beauty technology based on deep tissue stimulation. The gentle thermal action of RF (radiofrequency) denatures collagen fibers creating a lift effect and boosts the production of new collagen through the natural skin reparative mechanisms.

#### Technical Specifications

Voltage 100-230 Vac 50/60 Hz. Power Consumption Max 50W  
Power emission Max: 25W  
Frequency: 0,5MHz (body); 1MHz (face)  
Fuses: 2 x T2A (5 x 20 mm).  
Size: 180 x 380 x 380 mm Weight: 5 kg  
Handles: a tripolar handle and a hexapolar handle

### Multipolar Radiofrequency because...

The Monopole radiofrequency enters deep into the skin but it uses high power values and produces slow results. The Bipolar radiofrequency works on the surface, there is a high risk of burnings and it requires medium-high power value to get results. The **Tripolar or multipolar RF** is based on the fact that three or more electrodes direct the energy towards the centre among the three poles: this method is extremely safe and allows to address all the energy available in one point achieving a deep penetration at low power values (max 25-30 Watt) . The tripolar and hexapolar treatment is **pleasant like a warm massage** and it let you **work on skin laxity, wrinkles, local adiposities and cellulite**.

### Comfort and Effectiveness

**Firming Face Tripolar handle** studied for anti-wrinkle treatments and facial laxity treatments, it is perfect to work on the face, from the chin to the forehead, and to tone up neck and arms skin.

**Reshaping Body Hexapolar handle** designed to work on belly, hips, calves and thighs. To tone up, firm, fight skin laxity, reshape and let the client loose important cm at the thighs.



Powered by oxygen since 1990

Via B. Buozzi 8, 20060 Liscate (MI), Italy  
Tel 02 95350123/4 Fax 02 95350099  
[www.oxygensystem.net](http://www.oxygensystem.net) [www.o2bar.net](http://www.o2bar.net)  
[exar@o2bar.net](mailto:exar@o2bar.net)